


# Compatibility

# MICRODENT®

## HA Serie MICRODENT® SYSTEM

[ $\phi 3.5$  /  $\phi 4.2$  /  $\phi 5.1$ ]



### INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

#### Platform



$\phi 5.1$

IPD/NA-SW-00

IPD/NA-SW-01

Captive screw included. 1.20 Hex screwdriver.



### Ti-BASE Custom Ti-Base System®

#### Platform

#### Height (mm)



$\phi 5.1$

1.3

IPD/NA-IW-02/3D

IPD/NA-IW-03/3D



Please see [page 147](#) for the  Custom Ti-Base System® accessories.

For the Ti-Base it is recommended to use the screw IPD/NA-TW-00 or the IPD/MB-TR-50 for angled solutions.



### SCREW / TPA SCREW FOR ANGLED PROSTHESIS

| Platform              | Screwdriver | Seat    | Reference    | TiN Reference    |
|-----------------------|-------------|---------|--------------|------------------|
| $\phi 3.5$            | 1.20 Hex    | Conical | IPD/HA-TN-00 | IPD/HA-TN-00/TIN |
| $\phi 4.2$ $\phi 5.1$ | 1.75 Hex    | Conical | IPD/HA-TR-00 | IPD/HA-TR-00/TIN |
| $\phi 5.1$            | 1.25 Hex    | Conical | IPD/NA-TW-00 | -                |
| $\phi 3.5$            | TPA         | Conical | -            | IPD/HA-TN-50     |
| $\phi 4.2$ $\phi 5.1$ | TPA         | Conical | -            | IPD/HA-TR-50     |
| $\phi 5.1$            | TPA         | Conical | -            | IPD/MB-TR-50     |



### ANALOG

| Platform   | Reference    | Digital Reference |
|------------|--------------|-------------------|
| $\phi 3.5$ | IPD/HA-AN-00 | -                 |
| $\phi 4.2$ | IPD/HA-AR-00 | -                 |
| $\phi 5.1$ | IPD/HA-AW-00 | IPD/NA-AW-00/3D   |



TEMPORARY ABUTMENT / OPEN TRAY COPING

| Platform |  |  | Coping screw   |
|----------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|----------------|
| ø5.1     | IPD/NA-PW-06                                                                      | IPD/NA-PW-07                                                                      | IPD/NA-TW-00/C |

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

| Platform | Height (mm) | Reference    |
|----------|-------------|--------------|
| ø5.1     | 3.0         | IPD/NA-DW-03 |
|          | 5.0         | IPD/NA-DW-05 |



CASTABLE

| Platform |  |  |
|----------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
| ø3.5     | IPD/HA-HN-00                                                                      | IPD/HA-RN-00                                                                      |
| ø4.2     | IPD/HA-HR-00                                                                      | IPD/HA-RR-00                                                                      |
| ø5.1     | IPD/HA-HW-00                                                                      | IPD/HA-RW-00                                                                      |

