



Compatibility



SWEDEN & MARTINA®

MB Serie PREMIUM™ KOHNO®

[ϕ 3.3 / ϕ 3.8 / ϕ 4.25 / ϕ 5.0]





INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Platform	Height (mm)		
ϕ 3.3	10	IPD/MB-ST-00	IPD/MB-ST-01
	15	IPD/MB-ST-02	-
ϕ 3.8	10	IPD/MB-SN-00	IPD/MB-SN-01
	15	IPD/MB-SN-02	-
ϕ 4.25	10	IPD/MB-SR-00	IPD/MB-SR-01
	15	IPD/MB-SR-02	-

Captive screw included. 1.20 Hex Screwdriver.



Ti-BASE Custom Ti-Base System®


Platform	Height (mm)		
ϕ 3.3	0.5	IPD/MB-IT-00/3D	IPD/MB-IT-01/3D
	1.5	IPD/MB-IT-02/3D	IPD/MB-IT-03/3D
ϕ 3.8	0.7	IPD/MB-IN-00/3D	IPD/MB-IN-01/3D
	1.5	IPD/MB-IN-02/3D	IPD/MB-IN-03/3D
ϕ 4.25	0.7	IPD/MB-IR-00/3D	IPD/MB-IR-01/3D
	1.5	IPD/MB-IR-02/3D	IPD/MB-IR-03/3D

High version of the non-engaging Custom Ti-Base subject to availability.
Please see [page 147](#) for the  Custom Ti-Base System® accessories.



SCREW / TPA SCREW FOR ANGLED PROSTHESIS

Platform	Screwdriver	Seat	Reference	TiN Reference
ϕ 3.3 ϕ 3.8	1.25 Hex	Conical	IPD/MB-TN-00	IPD/MB-TN-00/TIN
ϕ 4.25 ϕ 5.0	1.25 Hex	Conical	IPD/MB-TR-00	IPD/MB-TR-00/TIN
ϕ 3.3 ϕ 3.8	TPA	Conical	-	IPD/MB-TN-50
ϕ 4.25 ϕ 5.0	TPA	Conical	-	IPD/MB-TR-50



ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ϕ 3.3	IPD/MB-AT-00	IPD/MB-AT-00/3D
ϕ 3.8	IPD/MB-AN-00	IPD/MB-AN-00/3D
ϕ 4.25	IPD/MB-AR-00	IPD/MB-AR-00/3D
ϕ 5.0	IPD/MB-AW-00	-
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01





STRAIGHT MULTI-UNIT ABUTMENT

Platform	Height (mm)	Reference
ø3.8	1.5	IPD/MB-MN-01
	3.0	IPD/MB-MN-03
	4.0	IPD/MB-MN-04

Our straight Multi-units can be used for single and multiple restorations.
Please see [page 157](#) for the Multi-unit system accessories.





TEMPORARY ABUTMENT / OPEN TRAY COPING

Platform			Coping screw
ø3.3	IPD/MB-PT-06	IPD/MB-PT-07	IPD/MB-TN-00/C
ø3.8	IPD/MB-PN-06	IPD/MB-PN-07	IPD/MB-TN-00/C
ø4.25	IPD/MB-PR-06	IPD/MB-PR-07	IPD/MB-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



CASTABLE

Platform		
ø3.8	IPD/MB-HN-00	IPD/MB-RN-00
ø4.25	IPD/MB-HR-00	IPD/MB-RR-00
ø5.0	IPD/MB-HW-00	IPD/MB-RW-00



PSD OVERDENTURE ABUTMENT

Platform	Height (mm)	Reference
ø3.3	1.0	IPD/MB-LT-01
	3.0	IPD/MB-LT-03
ø3.8	1.0	IPD/MB-LN-01
	3.0	IPD/MB-LN-03
	5.0	IPD/MB-LN-05
ø4.25	1.0	IPD/MB-LR-01
	3.0	IPD/MB-LR-03
	5.0	IPD/MB-LR-05



Please see [page 169](#) for the PSD system accessories.



HEALING ABUTMENT

Platform	Height (mm)	Reference
ø3.8	2.0	IPD/MB-DN-02
	3.0	IPD/MB-DN-03
	4.0	IPD/MB-DN-04
	6.0	IPD/MB-DN-06
	2.0	IPD/MB-DR-02
ø4.25	3.0	IPD/MB-DR-03
	4.0	IPD/MB-DR-04
	6.0	IPD/MB-DR-06



STRAIGHT AND ANGLED CASTABLE ABUTMENT WITH A Co-Cr BASE

Straight, 15° & 25°

Platform	Angulation		
ø3.3	Straight	IPD/MB-BT-00	IPD/MB-BT-01
	15°	IPD/MB-BT-10	IPD/MB-BT-11
	25°	IPD/MB-BT-20	IPD/MB-BT-21
ø3.8	Straight	IPD/MB-BN-00	IPD/MB-BN-01
	15°	IPD/MB-BN-10	IPD/MB-BN-11
	25°	IPD/MB-BN-20	IPD/MB-BN-21
ø4.25	Straight	IPD/MB-BR-00	IPD/MB-BR-01
	15°	IPD/MB-BR-10	IPD/MB-BR-11
	25°	IPD/MB-BR-20	IPD/MB-BR-21

Co-Cr Base for digital workflow	Scan Body
ø3.3 Straight / Angled	IPD/MB-BT-00/3D IPD/MB-BT-01/3D IPD/SC-AN-B1
ø3.8 Straight / Angled	IPD/MB-BN-00/3D IPD/MB-BN-01/3D IPD/SC-AN-B1
ø4.25 Straight / Angled	IPD/MB-BR-00/3D IPD/MB-BR-01/3D IPD/SC-AN-B1

All angled castable abutments include a TPA screw.

